

Crizal® Availability Chart

FOR EYECARE PROFESSIONALS

Monster scratches appear
on ordinary lenses

Not on Crizal Avancé lenses
with Scotchgard Protector

New Crizal Avancé™ lenses
with Scotchgard™ Protector

The most significant breakthrough
in scratch resistance ever

Crizal with **SCOTCHGARD**
AVANCÉ™ PROTECTOR
3M

Live Life in the Clear.™



CRIZAL® AVAILABILITY CHART

for Eyecare Professionals

Effective
FEB 10, 2009

2009 : Release.01

	PLASTIC 1.50				Airwear®				Thin&Lite® 1.60-1.74					Trivex®	Trivex® Transitions® VI	
	Clear	Polarized	Transitions® VI Gray	Transitions VI Brown	Clear	Polarized	Transitions VI Gray	Transitions VI Brown	1.60 Clear	1.60 Transitions VI	1.60 Polarized	1.67 Clear	1.67 Transitions VI			Thin&Lite 1.74 (**)
TINTABLE	✓				✓				✓			✓				
SINGLE VISION																
SV	•	•	•	•	•	•	•		•	•	•	•		•		
SFSV AS					•	•	•	•	•	•		•	•			
†Essilor Sport Solutions™					•	•										
2.0 Safety SV					•											
BIFOCALS																
FT28	•	•	•	•	•							•				
FT35	•		•	•												
FT45	•															
ROUND 22	•															
TRIFOCALS																
7x28	•		•	•												
8x35	•		•	•												
MID DISTANCE																
Interview™	•															
Nikon® Online™	•											•				
Prio® Computer Lens	•															
PROGRESSIVE LENSES																
VARILUX	Varilux® Physio®	•	•	•	•	•	•	•	•	•	•	•	•	•		
	**Varilux Physio 360°®	•	•	•	•	•	•	•	•	•	•	•	•	•		
	Varilux® Panamic®	•	•	•	•	•	•	•	•	•		•	•	•		
	Varilux Comfort®	•	•	•	•	•	•	•	•	•	•	•	•	•		
	**Varilux Comfort 360°™	•	•	•	•	•	•	•	•	•	•	•	•	•		
	Varilux® Ellipse®	•		•		•	•	•	•	•	•	•	•			
	**Varilux Ellipse 360°™	•		•		•	•	•	•	•	•	•	•	•		
	Varilux Liberty®	•		•		•	•									
DEFINITY	DEFINITY®	•‡	•	•	•	•‡	•	•	•	•		•	•		•	•
	DEFINITY SHORT™	•	•	•	•	•	•	•	•	•		•	•		•	•
ESSILOR PALS	Accolade®	•		•	•	•	•	•	•	•	•	•	•	•		
	**Accolade Freedom™	•		•	•	•	•	•	•	•	•	•	•			
	SmallFit®	•		•	•	•	•	•	•	•		•	•			
	Ovation®	•		•	•	•	•	•				•	•			
	Essilor Natural®	•		•		•	•		•							
	Adaptar®	•		•		•										
Nikon	Presio i 13™	•		•								•				
	Presio i 15™	•				•						•				

NEW

**All lenses with Thin&Lite 1.74 material will systematically include choice of Crizal® Alizé® or Crizal Avancé™ lenses with Scotchgard™ Protector.
 Varilux Physio 360°, Varilux Ellipse 360°, Varilux Comfort 360°, and Accolade Freedom Lenses include Crizal® anti-reflective of choice.
 † Essilor Sport Solutions has six different lens colors available, there is a lens to complement the unique lighting condition of most outdoor sporting activities.
 ‡ Only DEFINITY 1.50 and Airwear can be tinted.

For more information, please visit www.crizal.com



CRIZAL® AVAILABILITY CHART

for Eyecare Professionals

Effective
FEB 10, 2009

2009 : Release.01

		PLASTIC 1.50				MID	POLY			HI-INDEX					Trivex®	Younger Drivewear™ 1.50		
		Clear	Polarized	Transitions® VI Gray	Transitions VI Brown	1.56 Polarized	Clear	Transitions VI	Polarized	1.60 Clear	1.60 Transitions VI	1.60 Polarized	1.67 Clear	1.67 Transitions VI	1.67 Polarized			
TINTABLE		✓																
NEW	KBco.	SFSV		●				●							●			
		FT28		●														
		7x28		●														
NEW		EOS™ PAL		●												●		
		Fusion I PAL		●														
		Fusion II PAL							●									
	Specialty Lens Corporation	Specialty Select™ SV Spherical														●		
		Specialty Select™ SV Aspheric															●	
		Specialty Select™ FT28															●	
		Specialty Select™ short corridor															●	
		Specialty Select™ progressive															●	
		iRx Polaroid® SV		●			●		●			●				●		
		iRx Polaroid® FT28		●			●											
		iRx Polaroid® FT35		●			●											
		iRx Polaroid® 7x28		●			●											
	iRx Polaroid® 8x35		●			●												
	iRx Polaroid® SHORT					●		●										
	iRx Polaroid® PRO		●			●		●			●				●			
	Shamir	Genesis™	●		●	●				●	●		●	●		●		
		Genesis w/ Attitude™								●	●							
		Office™	●					●										
		Creation™	●		●			●		●	●		●	●				
NEW		Piccolo®	●		●			●		●	●		●	●				
	Sola	Compact Ultra™	●		●	●		●					●	●				
		XL®	●															
		VIP®	●		●	●												
		SOLA MAX®	●		●													
		Sola One®	●		●	●							●	●				
		Percepta™	●															
	AO	Compact™	●		●				●									
	Younger	SV Spherical														●	●	
		SV Aspheric														●		
		Image®	●	●	●	●										●	●	

For more information, please visit www.crizal.com



CRIZAL® AVAILABILITY CHART

for Eyecare Professionals

Effective
FEB 10, 2009

2009 : Release.01

		PLASTIC 1.50				MID	POLY			HI-INDEX							
		Clear	Polarized	Transitions® VI Gray	Transitions VI Brown	1.56 Polarized	Clear	Transitions VI	Polarized	1.60 Clear	1.60 Transitions VI	1.60 Polarized	1.67 Clear	1.67 Transitions VI	1.67 Polarized	Trivex®	Younger Drivewear™ 1.50
TINTABLE		✓															
X-Cel	SV																•
Signet Arm.	Kodak®	•															
	Kodak® Concise™	•					•						•				•
	Kodak® Precise™	•					•						•				•
Optima	Optima Resolution™						•										
	Optima Response™						•										
Zeiss	GT2	•		•	•			•					•	•			
	Gradal Top	•		•	•												
Seiko	SFSV												•				
	Proceed II												•				
	Proceed III												•				

SPECIALTY ITEMS: on UNCOATED LENSES ONLY

Slaboff, High Adds, Bi-Concave, etc...	•																
--	---	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

For more information, please visit www.crizal.com



FINISHED SINGLE VISION

New and Improved Crizal Avancé™ Lenses with Scotchgard™ Protector

Material	Power Range	Size
Orma® Flat	+4.00 to -6.00 up to -2.00 cyl.	70mm
	-2.00 cyl. available to a -4.00	
Airwear® (Aspheric) & (Spherical)	+4.00 to -8.00 up to -2.00 cyl.	65/70mm
	-2.00 cyl. available to a -6.00	
Airwear® Transitions® VI Gray (Aspheric) & (Spherical)	+4.00 to -6.00 up to -2.00 cyl.	65/70mm
	-2.00 cyl. available to a -4.00	
Thin&Lite® 1.60 (Aspheric)	+6.00 to -10.00 up to -2.00 cyl.	65/70/75mm
	-2.00 cyl. available to a -10.00	
Thin&Lite® 1.60 Transitions® VI Gray (Aspheric)	+4.00 to -6.00 up to -2.00 cyl.	65/70/75mm
	-2.00 cyl. available to a -6.00	
Thin&Lite 1.67 (Aspheric)	+6.00 to -10.00 up to -2.00 cyl.	65/70/75/80 mm
	-2.00 cyl. available to a -10.00*	
Thin&Lite 1.74 (Aspheric)	-3.00 to -12.00 up to -2.00 cyl.	70/75/80 mm
	-2.00 cyl. available to a -10.00	
1.50 Transitions VI Gray	+4.00 to -4.00 up to -2.00 cyl.	65/70 mm
	-2.00 cyl. available to a -2.00	

COMING
SOON

Crizal® Alizé®

Material	Power Range	Size
Airwear (Aspheric) & (Spherical)	+4.00 to -8.00 up to -2.00 cyl.	65/70mm
	-2.00 cyl. available to a -6.00	
Airwear Transitions VI Gray	+4.00 to -8.00 up to -2.00 cyl.	65/70mm
	-2.00 cyl. available to a -6.00	
Orma (Flat)	+4.00 to -6.00 up to -2.00 cyl.	70mm
	-2.00 cyl. available to a -4.00	
Thin&Lite 1.60 (Spherical)	-2.00 to -10.00 up to -2.00 cyl.	65/70/75mm
	-2.00 cyl. available to a -10.00	
Thin&Lite 1.60 (Aspheric)	-2.00 to -10.00 up to -2.00 cyl.	65/70/75mm
	-2.00 cyl. available to a -10.00	
Thin&Lite 1.67 (Spherical)	-2.00 to -8.00 up to -2.00 cyl.	70/75 mm
	-2.00 cyl. available to a -6.00	
Thin&Lite 1.67 (Aspheric)	-2.00 to -10.00 up to -2.00 cyl.	70/75/80 mm
	-2.00 cyl. available to a -10.00	

NEW

*Also available from -10.00 to -12.00 in Spherical only

For more information, please visit www.crizal.com



CRIZAL® STOCK AVAILABILITY CHART

for Eyecare Professionals

Effective
FEB 10, 2009

2009 : Release.01

FINISHED SINGLE VISION

Crizal®		
Material	Power Range	Size
Airwear® (Aspheric)	+4.00 to -8.00 up to -2.00 cyl.	65/70mm
	-2.00 cyl. available to a -6.00	
Airwear	+4.00 to -8.00 up to -2.00 cyl.	65/70mm
	-2.00 cyl. available to a -6.00	
Orma®	+4.00 to -6.00 up to -2.00 cyl.	70/75mm
	-2.00 cyl. available to a -4.00	
Orma Thin	+4.00 to -5.00 up to -2.00 cyl.	70mm
	-2.00 cyl. available to a -3.00	
Thin&Lite® 1.67 (Spherical)	-2.00 to -8.00 up to -2.00 cyl.	70/75mm
	-2.00 cyl. available to a -6.00	
Thin&Lite 1.67 (Aspheric)	-2.00 to -8.00 up to -2.00 cyl.	75/80mm
	-2.00 cyl. available to a -6.00	

For more information, please visit www.crizal.com

